

dpl® Gum Device User Guide updates needed to charging & battery.

CHARGING INSTRUCTIONS

Only use provided Power Supply. Attempts to use a third party Power Supply could damage device and will void the warranty.

- Charging:
- Charging typically takes 4-6 hours but may take longer.
- Separate the Oral Care Control Device from the Interchangeable Mouth Tray.
- Plug Power Cable into Power Supply and Oral Care Control Device.
- Plug Power Supply into power source.
- While charging, red light flashes on Oral Care Control Device on/off button.
- When fully charged, red light is solid red on the Oral Care Control Device on/off button.

BATTERY INDICATOR – Oral Care Control Device on/off button

Key to using the device is to make sure the LED's on the mouth guard tray light up.

- Battery Indicator:
- Check Oral Care Control Device on/off button.
- If **green** – device is charged and ready for use.
- If **red** – device needs to be charged.
- LED lights on the Mouth Guard Tray light up when on/off button on the Oral Care Control Device is pressed.
- To turn on device press and hold the on/off button for 2 seconds or until LED lights on the Mouth Guard Tray light up.
- dpl® Oral Care Gum Device will auto shut off after 10 minutes.
 - a. This time may vary based on the battery level of the Oral Care Control Device.